

Panther SL 637



APL Control LSO

The Panther SL spinal brace is indicated for the patient with moderate to severe low-back pain. The breathable mesh and the rigid Kydex posterior and lateral plates combine to provide bi-planar LSO control.

The string-and-pulley compression system allows the brace to be easily snugged, even by elderly arthritic hands.

The Panther SL is universal. It requires no tools to size. There's no need to cut or fold a belt. Just add or remove belt segments.

PDAC Approved: L0637/L0650



SCAN FOR
PANTHER CATALOG



The Best Value in Spinal Bracing

800-867-6070 | www.bluedortho.com

Panther SL 637



SPECIFICATIONS

- Bilateral pull handles
- Breathable-mesh belt
- Superior sagittal control
- Full-height posterior plate
- Conforms to any body shape
- Smooth dual-pull compression
- Lateral plates for coronal control

UNIVERSAL SIZE

- One size fits all.
- 25" to 66" and beyond! Just add or remove belt segments.

PDAC APPROVED

- L0637 Custom Fit
- L0650 Off the Shelf (OTS)

LORDOSIS

The posterior polymer plate is preset at 15°. It may be heat molded for exact fit during custom fitting.

